

INTUITIVE



# Da Vinci X & Da Vinci Xi Instrument & Accessory Catalogue

September 2022

Some products, features or technologies may not be available in all countries.  
Please contact your local Intuitive representative for product availability in your region.

For more information, call Intuitive.

+800 0821 2020 (toll-free) / +41 21 821 2020  
or visit [www.intuitive.com/en-gb](http://www.intuitive.com/en-gb)

For the name of your local distributor,  
call Intuitive or visit [www.intuitive.com/en-gb](http://www.intuitive.com/en-gb).

## Table of contents

### 8 mm EndoWrist instruments

- 4 EndoWrist monopolar cautery instruments
- 4 EndoWrist bipolar cautery instruments
- 6 EndoWrist clip appliers
- 7 EndoWrist needle drivers
- 8 EndoWrist graspers
- 9 EndoWrist scissors
- 10 Specialty EndoWrist instruments
- 10 Ultrasonic energy instrument
- 11 Vessel sealer extend
- 11 SynchroSeal
- 11 Suction irrigation

### Stapling components

- 12 SureForm stapler
- 12 SureForm stapler reloads
- 13 EndoWrist stapler 30
- 13 EndoWrist stapler 30 reloads
- 14 12 mm stapling accessories
- 14 Stapler accessories

### Accessories

- 15 Da Vinci X drapes
- 15 Da Vinci Xi drapes
- 15 Hasson cone
- 15 Reusable accessories
- 16 8 mm cannulas, obturators, and seals
- 16 8 mm disposable instrument accessories
- 16 Specialty accessories

### Vision equipment




- 17 8 mm endoscopes
- 17 Sterilization tray and connector set
- 17 Handheld camera components

### Energy accessories



- 18 Energy activation cables
- 18 Energy instrument cords
- 18 Da Vinci energy
- 18 Da Vinci simulation
- 18 Da Vinci Xi integrated table motion
- 19 Training equipment
- 21 Recommended spares and emergency equipment
- 21 System documentation
- 22 Technical support and customer service

## 8 mm instruments

### EndoWrist monopolar cautery instruments

	Part number
	<b>470179</b>
10 uses	
Shown with <b>tip cover accessory (1 procedure)</b>	<b>400180</b>
	<b>470184</b>
10 uses	
	<b>470183</b>
10 uses	

### EndoWrist bipolar cautery instruments

	Part number
	<b>471172</b>
14 uses	
	<b>471344</b>
14 uses	

8 mm  
instruments

EndoWrist bipolar cautery instruments (continued)

Part number



**Long bipolar grasper**  
14 uses

**471400**



**Fenestrated bipolar forceps**  
14 uses

**471205**



**Micro bipolar forceps**  
14 uses

**471171**






**Force bipolar**  
12 uses

**471405**

# 8 mm instruments

## EndoWrist clip appliers

		Part number
	<p><b>Medium-large clip applier<sup>1</sup></b>                      100 closures                      (Shown with Weck Hem-o-lok medium-large polymer clip)</p> <p>Validated for use only with Weck Hem-o-lok medium-large polymer clip.</p> <p><b>Weck product code:</b></p>	<p><b>470327</b></p> <p>544230</p>
	<p><b>Large clip applier<sup>2</sup></b>                      100 closures                      (Shown with Weck Hem-o-lok large polymer clip)</p> <p>Validated for use only with Weck Hem-o-lok large polymer clip.</p> <p><b>Weck product code:</b></p>	<p><b>470230</b></p> <p>544240</p>
	<p><b>Small clip applier<sup>3</sup></b>                      100 closures                      (Shown with Weck Horizon small-wide titanium clip)</p> <p>Validated for use only with Weck Horizon small-wide titanium clip.</p> <p><b>Weck product codes:</b></p>	<p><b>470401</b></p> <p>001201 001205</p>

1 Weck Hem-o-lok medium-large polymer clip (Weck product code 544230) is the only compatible clip used with this instrument. To order, contact Teleflex Medical ([www.teleflex.com](http://www.teleflex.com)) or your local distributor.

2 Weck Hem-o-lok large polymer clip (Weck product code 544240) is the only compatible clip used with this instrument. To order, contact Teleflex Medical ([www.teleflex.com](http://www.teleflex.com)) or your local distributor.

3 Weck Horizon small-wide titanium clip (Weck product code 001201 or 001205) is the only compatible clip used with this instrument. To order, contact TeleflexMedical ([www.teleflex.com](http://www.teleflex.com)) or your local distributor.

8 mm  
instruments

EndoWrist needle drivers

Part number



**Large needle driver**  
15 uses

**471006**



**Mega needle driver**  
10 uses

**470194**



**Mega SutureCut needle driver**  
15 uses

**471309**



**Large SutureCut needle driver**  
15 uses

**471296**

8 mm  
instruments


EndoWrist graspers

Part number



**ProGrasp forceps**  
18 uses

**471093**



**Long tip forceps**  
18 uses

**471048**



**Small grasping retractor  
(small graptor)**  
10 uses

**470318**



**Cobra grasper**  
18 uses

**471190**



**Tenaculum forceps**  
10 uses

**470207**



8 mm  
instruments

EndoWrist graspers (continued)

Part number



**Tip-Up fenestrated grasper**  
10 uses

**470347**



**Cadiere forceps**  
18 uses

**471049**

EndoWrist scissors

Part number



**Potts scissors**  
10 uses

**470001**



**Round tip scissors**  
10 uses

**470007**

## 8 mm instruments

### Specialty EndoWrist instruments

Part number



**Resano forceps**  
10 uses

**470181**



**Dual blade retractor**  
10 uses

**470249**



**Black diamond micro forceps**  
15 uses

**470033**



**DeBakey forceps**  
10 uses

**470036**

### Ultrasonic energy instruments

Part number



**Harmonic ACE curved shears**  
Box of 6  
single use

**480275**

8 mm  
instruments

Vessel sealer extend

Part number



**Vessel sealer extend**

Box of 6  
single use

**480422**

SynchroSeal

Part number



**SynchroSeal**

Box of 6  
single use

**480440**

Suction irrigation

Part number






**Suction irrigator**

Box of 6  
single use










**480299**

# Stapling components

## SureForm stapler



	Product description	Qty/ box	Uses	Part number
	SureForm 60 stapler	6	Single	<b>480460</b>
	SureForm 45 stapler	6	Single	<b>480445</b>
	SureForm 45 curved-tip stapler	6	Single	<b>480545</b>

## SureForm stapler reloads





	Product description	Qty/ box	Uses	Part number
	SureForm 60 white reload (2.5 mm, 6-row)	12	1 firing	<b>48360W</b>
	SureForm 60 blue reload (3.5 mm, 6-row)	12	1 firing	<b>48360B</b>
	SureForm 60 green reload (4.3 mm, 6-row)	12	1 firing	<b>48360G</b>
	SureForm 60 black reload (4.6 mm, 6-row)	12	1 firing	<b>48360T</b>
	SureForm 45 grey reload (2.0 mm, 6-row)	12	1 firing	<b>48345M</b>
	SureForm 45 white reload (2.5 mm, 6-row)	12	1 firing	<b>48345W</b>
	SureForm 45 blue reload (3.5 mm, 6-row)	12	1 firing	<b>48345B</b>
	SureForm 45 green reload (4.3 mm, 6-row)	12	1 firing	<b>48345G</b>
	SureForm 45 black reload (4.6 mm, 6-row)	12	1 firing	<b>48345T</b>

# Stapling components

## EndoWrist stapler 30









	Product description	Qty/ box	Uses	Part number
	EndoWrist stapler 30	1	50 firing	<b>470430</b>
	EndoWrist stapler 30 curved-tip	1	50 firing	<b>470530</b>

## EndoWrist stapler 30 reloads



	Product description	Qty/ box	Uses	Part number
	Stapler 30 grey reload (2.0 mm, 6-row)	12	1 firing	<b>48630M</b>
	Stapler 30 white reload (2.5 mm, 6-row)	12	1 firing	<b>48630W</b>
	Stapler 30 blue reload (3.5 mm, 6-row)	12	1 firing	<b>48630B</b>
	Stapler 30 green reload (4.3 mm, 6-row)	12	1 firing	<b>48630G</b>

## Stapling components

### 12 mm stapling accessories

	Product description	Qty/box	Uses	Part number
	12 mm and stapler cannula (100 mm)	1	NA	<b>470375</b>
	12 mm and stapler blunt obturator	1	NA	<b>470376</b>
	12 mm and stapler cannula, long (150 mm)	1	NA	<b>470389</b>
	12 mm and stapler blunt obturator, long	1	NA	<b>470390</b>
	12 mm and stapler cannula seal (box of 10)	10	1 procedure	<b>470380</b>
	12 - 8 mm reducer (box of 6)	6	1 procedure	<b>470381</b>
	12 mm and stapler bladeless obturator	1	NA	<b>470395</b>
	12 mm and stapler bladeless obturator, long	1	NA	<b>470396</b>

### Stapler accessories

	Product description	Qty/box	Uses	Part number
	Stapler sheath (box of 10)	10	1 procedure	<b>410370</b>
	Stapler release kit <sup>1</sup>	1	NA	<b>381215</b>

<sup>1</sup> Documentation is listed in English only. For other languages contact customer service. Documentation may not be available in all languages.

## Accessories



## Da Vinci X drapes

Product description	Qty/ box	Uses	Part number
Arm drape (box of 20)	20	1 procedure	<b>470015</b>
Arm-4 extension drape (box of 20)	20	1 procedure	<b>470473</b>

## Da Vinci Xi drapes

Product description	Qty/ box	Uses	Part number
Arm drape (box of 20)	20	1 procedure	<b>470015</b>
Column drape (box of 20)	20	1 procedure	<b>470341</b>

## Hasson cones

Product description	Qty/ box	Uses	Part number
 8 mm hasson cone	1	NA	<b>470398</b>
 12 mm hasson cone	1	NA	<b>470399</b>









## Reusable accessories

Product description	Qty/ box	Uses	Part number
Instrument release kit (IRK) <sup>1</sup>	1	NA	<b>381216</b>
Gage pin (box of 3)	3	NA	<b>470397</b>
Blue fiber cable kit	1	NA	<b>380989</b>
Instrument introducer	1	NA	<b>342562</b>


<sup>1</sup> Documentation is listed in English only. For other languages contact customer service. Documentation may not be available in all languages.

## Accessories




## 8 mm cannulas, obturators, and seals

	Product description	Qty/ box	Uses	Part number
	5 mm - 8 mm cannula seal (box of 10)	10	1 procedure	<b>470361</b>
	8 mm cannula	1	NA	<b>470002</b>
	8 mm cannula, long	1	NA	<b>470004</b>
	8 mm flared/grounded cannula	1	NA	<b>470319</b>
	8 mm blunt obturator	1	NA	<b>470008</b>
	8 mm blunt obturator, long	1	NA	<b>470009</b>
	8 mm bladeless obturator (optical) (box of 6)	6	1 procedure	<b>470359</b>
	8 mm bladeless obturator, long (optical) (box of 6)	6	1 procedure	<b>470360</b>

## 8 mm disposable instrument accessories

	Product description	Qty/ box	Uses	Part number
	Tip cover accessory (box of 10)	10	1 procedure	<b>400180</b>



## Specialty accessories

	Product description	Qty/ box	Uses	Part number
	8 mm camera cannula	1	NA	<b>478263</b>
	8 mm blunt obturator	1	NA	<b>478008</b>
	5 mm-10 mm cannula seal (box of 10)	10	1 procedure	<b>478161</b>





# Vision equipment





## 8 mm endoscopes

	Product description	Qty/ box	Uses	Part number
	8 mm endoscope plus, 0°	1	NA	<b>470056</b>
	8 mm endoscope plus, 30°	1	NA	<b>470057</b>

## Sterilization tray and connector set



	Product description	Qty/ box	Uses	Part number
	8 mm endoscope sterilization tray	1	NA	<b>400498</b>
	8 mm endoscope washer/ disinfector connector set	1	NA	<b>400482</b>

## Handheld camera components



	Product description	Qty/ box	Uses	Part number
	Handheld camera	1	NA	<b>470035</b>
	Handheld camera light guide cable	1	NA	<b>470042</b>
	Handheld camera light guide adapter	1	NA	<b>470051</b>
	Handheld camera stainless steel reprocessing tray	1	NA	<b>400499</b>

## Energy accessories


### Energy activation cables

	Product description	Qty/ box	Uses	Part number
	Energy activation cable, Covidien, Force Triad	1	NA	<b>371716</b>
	Energy activation cable, Ethicon, Gen 11	1	NA	<b>371870</b>





### Energy instrument cords

	Product description	Qty/ box	Uses	Part number
	Monopolar energy instrument cord (4 m)	1	NA	<b>470383</b>
	Bipolar energy instrument cord (5 m)	1	NA	<b>470384</b>


## Da Vinci energy

	Product description	Qty/ box	Uses	Part number
	E-100 generator	1	NA	<b>378848</b>

## Da Vinci simulation

	Product description	Qty/ box	Uses	Part number
	Da Vinci SimNow simulator <sup>1</sup>	1	NA	<b>377773</b>
	Da Vinci SimNow service and benefits (annual)	1	NA	<b>600092</b>
	Blue fiber cable, simulator (1 m)	1	NA	<b>381129</b>
	Simulator skills drills	1	NA	<b>372363</b>

## Da Vinci Xi integrated table motion

	Product description	Qty/ box	Uses	Part number
	Integrated table motion hardware and TAM table access software upgrade <sup>2</sup>	1	NA	<b>600062</b>

<sup>1</sup> The da Vinci SimNow simulator requires an annual service fee (PN 600092). For details contact Intuitive customer service.

<sup>2</sup> The TS 7000dV table and dV remote control are sold and serviced by Hillrom. To order, contact Hillrom or your local distributor (<https://www.hill-rom.com/international/Our-Company/Contact-Us/Hill-Rom-Office-Locations/>).

# Training equipment

## Training instruments

Product description	Uses	Part number
Large needle driver	30	<b>470006-T</b>
ProGrasp forceps	30	<b>470093-T</b>
Tip-Up fenestrated grasper	30	<b>470347-T</b>
Long bipolar grasper	30	<b>470400-T</b>
Maryland bipolar forceps	30	<b>470172-T</b>
Monopolar curved scissors <sup>1,2</sup>	30	<b>470179-T</b>
Permanent cautery hook	30	<b>470183-T</b>
Permanent cautery spatula	30	<b>470184-T</b>
Fenestrated bipolar forceps	30	<b>470205-T</b>
Force bipolar	30	<b>470405-T</b>
Mega SutureCut needle driver	30	<b>470309-T</b>
Large clip applier <sup>3</sup>	200 closures	<b>470230-T</b>
Small grasping retractor (small graptor)	30	<b>470318-T</b>
Curved bipolar dissector	30	<b>470344-T</b>
Medium-large clip applier <sup>4</sup>	200 closures	<b>470327-T</b>
Cadiere forceps	30	<b>470049-T</b>
Large SutureCut needle driver	30	<b>470296-T</b>
Mega needle driver	30	<b>470194-T</b>
Small clip applier <sup>5</sup>	200 closures	<b>470401-T</b>

# Training equipment

## Training instruments (continued)

Product description	Uses	Part number
EndoWrist stapler 30	125 firings	<b>470430-T</b>
EndoWrist stapler 30 curved-tip	125 firings	<b>470530-T</b>
SureForm 60 stapler (box of 6)	24 firings	<b>480460-T</b>
SureForm 45 stapler (box of 6)	24 firings	<b>480445-T</b>
SureForm 45 curved-tip stapler (box of 6)	24 firings	<b>480545-T</b>
EndoWrist stapler 30 demo reload (Box of 20)	1	<b>382051-T</b>

1 Requires tip cover accessory (PN 400180) for use

2 Requires in-service prior to first use, please contact your local Intuitive Sales Representative.

3 Weck Hem-o-lok large polymer clip (Weck product code 544240) is the only compatible clip used with this instrument.

To order, contact Teleflex Medical ([www.teleflex.com](http://www.teleflex.com)) or your local distributor.

4 Weck Hem-o-lok medium-large polymer clip (Weck product code 544230) is the only compatible clip used with this instrument.

To order, contact Teleflex Medical ([www.teleflex.com](http://www.teleflex.com)) or your local distributor.

5 Weck Horizon small-wide titanium clip (Weck product code 001201 or 001205) is the only compatible clip used with this instrument.

To order, contact Teleflex Medical ([www.teleflex.com](http://www.teleflex.com)) or your local distributor.

Notes: (a) Training instruments are not approved for human use. (b) Training instruments can be readily identified by the red instrument housing. (c) Red instrument icons persist on the monitors when a training instrument is loaded onto the system. (d) The training instrument warranty covers only defects in manufacturing. Upon first use, the warranty expires.

## Recommended spares and emergency equipment

Product description	Page number	Part number
Endoscopes	17	Multiple
Blue fiber cable kit	15	380989
EndoWrist instruments	4-9	Multiple
Instrument cannula	14, 16	Multiple
Instrument release kit (IRK) (item provided free of charge)	15	381216
Gage pin (set of 3) (item provided free of charge)	15	470397

## System documentation

Product description	Part number
System documentation kit – includes: • User manuals • Quick reference guides (QRG)	381227*
Reprocessing documentation kit – includes: • Reprocessing instructions for instruments, accessories and endoscopes	470601*
Universal reprocessing hardware kit – includes: • Tip cover removal pictogram • Luer fitting for flushing • Pipe-to-hose adapter • Brush with soft, nylon bristles • Tip cover insertion/removal tool for monopolar curved scissors (Qty. 5)	470629

\*Documentation is listed in English only. For other languages contact Intuitive customer service. Documentation may not be available in all languages.

Note: Content and part numbers of the kits may vary. More documentation is available on demand and in PDF format. For details contact Intuitive customer service.

# Technical Support and Customer Service

## Contact us by phone

+800 0821 2020 (toll free) or use a local number:

<b>Germany</b> +49 698 991 4175	<b>Austria</b> +43 1 26 75 225	<b>Switzerland</b> +41 21 821 2020	<b>France</b> +33 1 77 68 88 45
<b>Belgium</b> +32 2 588 01 62	<b>The Netherlands</b> +31 1 08 08 03 28	<b>Sweden</b> +46 40 6450546	<b>Denmark</b> +45 78772514
<b>Norway</b> +47 21984147	<b>Finland</b> +358 9 31579155	<b>United Kingdom</b> +44 2 03 37 62 782	<b>Ireland</b> +353 1 69 17 094

## The options during your call

### Option 1: Technical Support

#### Emergency technical support

Technical issues with Intuitive products before or during surgery.

**24/7 available** (in local language\*)

Mon. to Fri. 7am to 5pm (CET),  
Afterhours support only in English.

\* English, French, German,  
Dutch and Spanish.

### Option 2: Technical Support

#### Technical questions

Questions related to Intuitive products and technology.

**Available Mon. to Fri.  
7am to 5pm (CET).**

### Option 3: Customer Service

#### Orders and returns

- Purchase order inquiries
- Order delivery status
- Complaint
- Product return authorization
- Endoscope, camera head & simulator exchanges

**Available Mon. to Fri.  
8am to 6pm (CET).**

You can access these services also 24/7 on the Intuitive Customer Portal ([www.davincisurgerycommunity.com](http://www.davincisurgerycommunity.com)). If you want a user account, please contact Intuitive customer service.

## Contact Customer Service by e-mail

### Germany:

Support.DE@intusurg.com

### Austria:

Support.AT@intusurg.com

### Switzerland:

Support.CH@intusurg.com

### France:

Support.FR@intusurg.com

### Belgium:

Support.BE@intusurg.com

### The Netherlands:

Support.NL@intusurg.com

### Luxembourg:

Support.LU@intusurg.com

### Sweden:

Support.SE@intusurg.com

### Denmark:

Support.DK@intusurg.com

### Norway:

Support.NO@intusurg.com

### Finland:

Support.FI@intusurg.com

### United Kingdom:

Support.UK@intusurg.com

### Ireland:

Support.IE@intusurg.com

### Indirect Markets:

Support.Indirect@intusurg.com

## European Regional Offices

### European Headquarters (Switzerland)

Intuitive Surgical Sàrl  
1, chemin des Müriers  
1170 Aubonne, Suisse  
Tel: +41 21 821 2000  
Fax: +41 21 821 2001  
VAT #: CHE-113.407.162 VAT

### Germany

Intuitive Surgical Deutschland GmbH  
Am Flughafen 6  
79108 Freiburg, DE  
MWST #: DE200038727  
Amtsgericht Freiburg, HRB 720585  
Geschäftsführer: Dirk Barten,  
Marshall Mohr,  
Jean-Yves Raimon-Dacunha-Castelle,  
Kara Andersen Reiter

### France

Intuitive Surgical SAS  
Cité de la Photonique,  
bâtiment Gienah  
11 avenue de Canteranne  
33600 Pessac, FR  
TVA #: FR 64 417 889 177  
SIREN 417 889 177  
Capital Social 37000€  
RCS Bordeaux Code APE 4646Z

### United Kingdom

Intuitive Surgical Ltd  
The Schrodinger Building  
Oxford Science Park  
Heatley Rd, Oxford, OX4 4GE, UK  
VAT #: GB943 2213 49  
Company registration  
No 6727762

## Global Headquarters

### USA

Intuitive Surgical, Inc.  
1020 Kifer Road, Building 108  
Sunnyvale, CA 94086-5304  
Tel: +1.408.523.2100  
Fax: +1.408.523.1390

For more information, call Intuitive.

+800 0821 2020 (toll-free) / +41 21 821 2020  
or visit [www.intuitive.com/en-gb](http://www.intuitive.com/en-gb)

For the name of your local distributor,  
call Intuitive or visit [www.intuitive.com/en-gb](http://www.intuitive.com/en-gb).

#### **Legal Notices**

Individuals' outcomes may depend on a number of factors, including but not limited to patient characteristics, disease characteristics and/or physician/surgeon experience.

The information contained in this document has been checked and compiled with the greatest care. However, no responsibility is taken for its correctness, completeness and topicality. It is the sole responsibility of the recipient to check all information before using it in the individual case.

Some products, features or technologies may not be available in all countries. Please contact your local Intuitive representative for product availability in your region.

Refer to the product specific User Manual for indications, contraindications, warnings and other product information.

#### **Privacy Notice**

Intuitive's Privacy Notice is available at [www.intuitive.com/privacy](http://www.intuitive.com/privacy).

© 2022 Intuitive Surgical Operations, Inc.  
All rights reserved. Product and brand names/ logos are trademarks or registered trademarks of Intuitive Surgical or their respective owner.